
# FINNPIPETTE COLOUR MULTICHANNEL FINNPIPETTE MULTICANAUX DIGITAL-COLORS

INSTRUCTIONS FOR USE  
GUIDE D'UTILISATION

# Thermo

SCIENTIFIC

Thermo Fisher Scientific Oy  
P.O.Box 100, Fi-01621 Vantaa, Finland  
Tel. +358-9-329 100, fax -358-9-3291 0414  
[www.thermo.com/finnpipette](http://www.thermo.com/finnpipette)

## ENGLISH

|                            |      |
|----------------------------|------|
| PRODUCT DESCRIPTION .....  | 3    |
| PIPETTE OPERATION .....    | 3    |
| PIPETTING TECHNIQUES ..... | 4    |
| MAINTENANCE .....          | 5    |
| CALIBRATION .....          | 6    |
| PACKAGE .....              | 8    |
| TROUBLE SHOOTING .....     | 8    |
| SPARE PARTS .....          | 9-10 |
| FIGURES .....              | 18   |

## FRANÇAIS

|                                 |       |
|---------------------------------|-------|
| DESCRIPTION .....               | 12    |
| UTILISATION DE LA PIPETTE ..... | 12    |
| METHODES DE PIPETAGE .....      | 13    |
| MAINTENANCE .....               | 14    |
| CALIBRAGE .....                 | 15    |
| CONDITIONNEMENT .....           | 17    |
| ERREURS ET CAUSES .....         | 17    |
| PIECES DETACHEES .....          | 10-11 |
| FIGURES .....                   | 18    |

## PRODUCT DESCRIPTION

FINNPIPETTE COLOUR MULTICHANNELS are handheld, adjustable volume, general purpose micropipettes for sampling and dispensing of accurate liquid volumes. They operate on the air displacement principle (i.e. an air interface) and use detachable, disposable tips or tip bands.

FINNPIPETTE COLOUR MULTICHANNELS are available as either 8- 12-channel models with a tip ejector system. Each model can be supplied in either of two volume ranges covering a total of 5 to 300  $\mu\text{l}$ . The two volume ranges are:

| Order No | Channel | Volume rang          | Scale | Increment         |
|----------|---------|----------------------|-------|-------------------|
| 4173200  | 8       | 5-50 $\mu\text{l}$   |       | 0,5 $\mu\text{l}$ |
| 4173210  | 8       | 50-300 $\mu\text{l}$ |       | 5,0 $\mu\text{l}$ |
| 4173220  | 12      | 5-50 $\mu\text{l}$   |       | 0,5 $\mu\text{l}$ |
| 4173230  | 12      | 50-300 $\mu\text{l}$ |       | 5,0 $\mu\text{l}$ |

## PIPETTE MATERIALS

FINNPIPETTE COLOUR MULTICHANNELS are made of mechanically durable and chemically resistant materials.

|                    |    |                   |    |
|--------------------|----|-------------------|----|
| Operating Cap      | PP | Tip ejector       | PA |
| Operating Button   | PA | Tip cones         | PC |
| Handle             | PA | Pipette body      | PA |
| Tip ejector pusher | PA | Tip ejector plate | PC |

## DESCRIPTION OF TIPS

Finntips are recommended for use with the FINNPIPETTE COLOUR MULTICHANNEL.

They are made of natural colour polypropylene, generally regarded as the only contamination free material suitable for tips. Finntips are also autoclavable (121°C).

Following accessories are useful with handling of pipettes and pipetting of reagents:

|              |         |
|--------------|---------|
| Stand        | 9420290 |
| Shelf hanger | 2206040 |

## DIGITAL DISPLAY (Fig. 1)

The adjusted delivery volume is clearly indicated by a digital display on the pipette handle.

## PIPETTE OPERATION

### SETTING THE DELIVERY VOLUME (Fig. 2)

The delivery volume is set using the operating button on the top of the pipette. To increase the delivery volume turn the operating button counterclockwise. To decrease turn it clockwise.

Ensure that the desired delivery volume clicks into place and that the digits for the volume are completely visible on the handle display.

Do not use force to try to set volumes outside of the pipette's specified volume range!

**TIP EJECTION (Fig. 3)**

Each pipette is fitted with a tip ejector system to help eliminate the possibility of contamination. For tip ejection, two methods are available:

- A.** If one-handed operation is preferred, twist the pipette handle so that the tip ejector pusher and tip ejector on the body assembly are in line. Press the tip ejector pusher with your thumb directing the tips towards a suitable waste receptacle.
- B.** In two-handed operation, press the tip ejector on the body assembly directing the tips towards a suitable waste receptacle. With some tip types this may be an easier and lighter method for user if the tips are very tight.

**PIPETTING TECHNIQUES**

**Figures 4-6:**

- A = Ready position**  
**B = First stop**  
**C = Second stop**

**FORWARD TECHNIQUE (Fig. 4)**

Ensure that clean tips are firmly pushed on to the tip cones of the pipette.

The operation of the pipette is controlled by the operating button. For the best possible accuracy:

Operate the button slowly and with steady speed at all times, particularly when working with liquids having high viscosity. Never let the operating button snap back. When working with serum or biological fluids prime the tips with the liquid before pipetting. This is done by filling and emptying the tips.

Hold the pipette vertically (10 degrees maximum from vertical) during the liquid intake. For maximum hand control, hold the pipette tightly in the hand with the "Grippy" resting on the index finger.

Fill a clean reagent reservoir with the liquid to be dispensed.

- 1.** Depress the operating button to the first stop.
- 2.** Dip the tips just under the surface of the liquid in the reservoir (2-3 mm) and slowly release the operating button. Withdraw the tips from the liquid touching them against the edge the reservoir to remove excess liquid.
- 3.** Deliver the liquid by gently depressing the operating button to the first stop. After a delay of about one second, continue to depress the operating button all the way down to the second stop. This action will empty the tips.
- 4.** Release the operating button to the ready position.

If necessary, change the tips and continue with the pipetting.

**REVERSE TECHNIQUE (Fig. 5)**

The reverse technique is suitable for dispensing liquids having high viscosity or a tendency to foam easily. The technique is also recommended for dispensing very small volumes. Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button all the way down to the second stop.
2. Dip the tips just under the surface of the liquid in the reservoir (2-3 mm) and slowly release the operating button.  
This action will fill the tips. Withdraw the tips from the liquid touching them against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the operating button to the first stop. Hold the operating button at the first stop. Some liquid will remain in the tip and should not be included in the delivery.
4. The remaining liquid is either discarded with the tips or pipetted back into the container.

### **REPETITIVE TECHNIQUE (Fig. 6)**

The repetitive technique offers a rapid and simple procedure for repeated delivery of the same volume. Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button all the way down to the second stop.
2. Dip the tips just under the surface of the liquid in the reservoir (2-3 mm) and slowly release the operating button.  
This action will fill the tips. Withdraw the tips from the liquid touching them against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the operating button to the first stop. Hold the operating button at the first stop. Some liquid will remain in the tip and should not be included in the delivery.
4. Continue pipetting by repeating the steps 2 and 3.

## **MAINTENANCE**

When the FINNPIPETTE COLOUR MULTICHANNEL is not in use make sure it is safely stored in a vertical position, preferably on a Finnpiquette Stand.

### **SHORT TERM CHECKING**

The pipette should be checked at the beginning of each day for dust and dirt on outside surfaces of the pipette. Particular attention should be paid to the tip cones. Solvents should not be used for cleaning the pipette.

### **LONG TERM MAINTENANCE**

If the pipette is used daily it should be checked at least twice a year. The servicing procedure is as follows:

1. Unscrew the two screws on the tip ejector plate and remove it.
2. Unscrew the eight screws on the body assembly cover and remove it.
3. Remove the tip ejector with the spring.
4. Depress the operating button all the way down. Remove the tip cone holder with the tip cones.

5. Clean the pistons and O-rings with a dry napless cloth and grease them slightly with grease (included in the package).
6. Reassemble the pipette.  
In the event of a malfunction clean the pipette and send it to the Thermo Fisher Scientific representative.  
Please, enclose a complete list of liquids which have entered the pipette.

## CALIBRATION

All Finnpiettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water. The pipettes are constructed to permit re-adjustment for liquids of different temperature and viscosity.

### DEVICE REQUIREMENTS AND TEST CONDITIONS

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

| Volume range            | readable graduation |
|-------------------------|---------------------|
| under 10 $\mu\text{l}$  | 0.00 1 mg           |
| 10-100 $\mu\text{l}$    | 0.01 mg             |
| above 100 $\mu\text{l}$ | 0.1 mg              |

Test liquid: Water, distilled or deionized, "grade 3" water conforming ISO 3696. Tests are done in a draft-free room at a constant ( $\pm 0.5^\circ\text{C}$ ) temperature of water, pipette and air between  $15^\circ\text{C}$  to  $30^\circ\text{C}$ .

The relative humidity must be above 50%. Especially with volumes under 50  $\mu\text{l}$  the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

#### Procedure to check calibration

The pipette is checked with the maximum volume (nominal volume) and with the minimum volume. A new tip is first pre-wetted 3-5 times and a series of ten pipettings is done with both volumes. A pipette is always adjusted for delivery (Ex) of the selected volume.

#### Procedure:

1. Do 10 pipettings with the minimum volume.
2. Do 10 pipettings with the maximum volume.
3. Calculate the inaccuracy (A) and imprecision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the calculated results are within the selected limits, the adjustment of the pipette is correct.

TABLE 1: Maximum permissible errors according ISO8655

**PACKAGE**

The FINNPIPETTE COLOUR MULTICHANNELS are shipped in a custom made, cardboard package.

The package contains the following items:

1. FINNPIPETTE MCP
2. Instruction Manual
3. Tube of Grease
4. Service tool
5. Samples of Tips

**TROUBLE SHOOTING**

The following table is a guide to possible problems and their solutions.


| <b>Fault</b>                               | <b>Possible reason</b>  | <b>Solution</b>  |
|--|---|--|
| Leakage                                    | Tip incorrectly attached  | Attach firmly  |
|  | Foreign bodies between tips and tip cones                                 | Clean tip cones thoroughly and attach new tips                           |
|  | Foreign bodies between O-rings and pistons                                | Clean and grease O-rings and the pistons in the cylinders.<br>Use grease |
|  | Insufficient amount of grease on pistons and O-rings                      | Grease accordingly   |
| Inaccurate dispensing                      | O-rings damaged   | Change the O-rings   |
|  | Incorrect operation   | Follow instructions carefully  |
|  | Tips incorrectly attached   | Attach firmly  |
| Inaccurate dispensing with certain liquids | Calibration altered: e.g. caused by misuse                                | Recalibrate according to instructions                                    |
|  | Unsuitable calibration. Many high viscosity liquids require recalibration | Recalibrate with the liquids in question                                 |

**CAUTION!**

The Finnpiquette is designed to allow easy in-lab service. If you would prefer to have us or your local representative service your pipette, please make sure that the pipette has been decontaminated before you send it to us. Please note that the postal authorities in your country may prohibit or restrict the shipment of contaminated material by mail.

| <b>Description</b>  |  |  | <b>Code</b>                                  |
|---|--|--|--|
| 1. Handle (incl. item 18)   |  | 5-50 µl<br>50-300 µl                         | 2205130<br>2205140                           |
| 2. Tip ejector  |  | 5-50 µl<br>50-300 µl                         | 10590140<br>10590150                         |
| 3. Tip ejector spring<br>4. Body assembly cover screw<br>5. Tip ejector screw |  |  | 1130111<br>0202061<br>0202021                |
| 6. Tip ejector plate  | 8-channel<br>8-channel<br>12-channel<br>12-channel | 5-50 µl<br>50-300 µl<br>5-50 µl<br>50-300 µl | 10590300<br>10590310<br>10590320<br>10590330 |
| 7. Body assembly cover  | 8-channel<br>8-channel<br>12-channel<br>12-channel | 5-50 µl<br>50-300 µl<br>5-50 µl<br>50-300 µl | 10590240<br>10590250<br>10590260<br>10590270 |
| 8. Tip cone   |  | 5-50 µl<br>50-300 µl                         | 10590340<br>10590350                         |
| 9. O-ring   |  | 5-50 µl<br>50-300 µl                         | 1030480<br>1030260                           |
| 10. Collar  | 8-12-channel<br>8-12-channel                       | 5-50 µl<br>50-300 µl                         | 1056030<br>1056040                           |
| 11. Tip cone holder   | 8-channel<br>12-channel<br>8-channel<br>12-channel | 5-50 µl<br>5-50 µl<br>50-300 µl<br>50-300 µl | 1056090<br>1056100<br>1056120<br>1056130     |
| 12. Label<br>13. Tip ejector pusher holder                                    |  |  | 1527110<br>1054820                           |
| 14. Spring support<br>15. Tip ejector pusher spring                           |  |  | 1056630<br>1130570                           |
| 16. Tip ejector pusher  |  | 5-50 µl<br>50-300 µl                         | 10588776<br>10588786                         |
| 17. Service tool<br>18. Cap   |  |  | 10589070<br>10588850                         |
| 19. Body  | 8-channel<br>8-channel<br>12-channel<br>12-channel | 5-50 µl<br>50-300 µl<br>5-50 µl<br>50-300 µl | 2205090<br>2205100<br>2205110<br>2205120     |
| 20. Piston holder   | 8-channel<br>12-channel                            |  | 10590770<br>10590810                         |
| 21. Holder spring   | 8-channel<br>12-channel                            |  | 1130680<br>1130720                           |
| 22. Locking pieces  |  |  | 1058180                                      |
| 23. Piston  |  | 5-50 µl<br>50-300 µl                         | 2201240<br>2202940                           |
| 24. Piston spring   | 8-12-channel                                       |  | 1131630                                      |

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Finnpipette® and Finntip® are registered trademarks of Thermo Scientific.
© 2007 Thermo Fisher Scientific Inc. All rights reserved.