# Dispensette® III Dispensette® Organic

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## Operating Manual

Please read carefully before use, and follow all operating and safety instructions!

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Read carefully before use!

## Safety Instructions

This instrument may be involved with hazardous materials, operations, and equipment.

This manual does not purport to address all of the safety risks associated with its use. It is the responsibility of whomever uses this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Follow the general instructions for hazard prevention and safety regulations, e.g., wear protective clothing, eye protection and gloves.
2. Observe the reagent manufacturer's information.
3. Every user must be acquainted with this Operating Manual before operation.
4. Use the instrument only for dispensing liquids, with strict regard to the defined Operating Exclusions and Limitations (► pages 7-8). If in doubt, contact the manufacturer or supplier.
5. When dispensing flammable media, make provisions to avoid electrostatic charging, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

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## Application and design

This instrument is designed for dispensing liquids, observing the following limits:

- 15 to 40 °C of instrument and reagent
- vapor pressure up to 500 mbar
- density up to 2.2 g/cm<sup>3</sup>
- viscosity up to 500 mm<sup>2</sup>/s  
(dynamic viscosity [mPas] = kinematic viscosity [mm<sup>2</sup>/s] x density [g/cm<sup>3</sup>])

## Materials

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:

Borosilicate glass, Al<sub>2</sub>O<sub>3</sub>, ETFE, PFA, FEP, PTFE,  
platinum-iridium (Dispensette® III)  
or tantalum (Dispensette® Organic),  
PP (screw cap).

## Selecting the right dispenser

See “Dispenser selection chart” (► page 41) for guidelines on selecting the right dispenser.

### Note:

Consult the following chapters for operating exclusions and operating limitations with your application (► pages 8, 9).

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### ■ Dispensette® III

#### Never use this instrument with

- liquids attacking  $\text{Al}_2\text{O}_3$ , ETFE, FEP, PFA and PTFE (e.g., dissolved sodium azide)
- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- liquids which are decomposed catalytically by platinum-iridium (e.g.,  $\text{H}_2\text{O}_2$ )
- fuming acids
- tetrahydrofuran
- trifluoroacetic acid
- Explosive liquids (e.g. carbon disulfide)
- suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument

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### ■ Dispensette® Organic

#### Never use this instrument with

- liquids attacking  $\text{Al}_2\text{O}_3$ , tantalum, ETFE, FEP, PFA and PTFE (e.g., dissolved sodium azide)
- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- bases and saline solutions
- Explosive liquids (e.g. carbon disulfide)
- suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument.

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## Operating Limitations

Liquids which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or highly concentrated alkaline solutions).

When dispensing flammable media, make provisions to avoid electrostatic charging, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

### Warning!

If there is a sign of a potential malfunction (e.g., piston difficult to move) never use force. Immediately stop dispensing and follow cleaning instructions (► page 19) or contact the manufacturer.

### Note:

Compatibility of the instrument for this special application (e.g., trace material analysis) must be checked by the user or contact the manufacturer.

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



► page 37.

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► pág. 37.

Model Digital • Easy Calibration is  
manufactured under U.S. Patent  
5,957,330.

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\*\* SafetyPrime™ Ventil (optional)  
erfordert Rückdosierrohr

## Components

|              |                                                         |
|--------------|---------------------------------------------------------|
| A            | Piston seat                                             |
| B            | Cylinder sleeve                                         |
| C            | Volume adjustment                                       |
| D            | Tension ring                                            |
| <b>E*</b>    | <b>Piston</b>                                           |
| <b>F*</b>    | <b>Glass cylinder</b>                                   |
| G            | Air vent cap                                            |
| <b>H*</b>    | <b>Valve block with serial number</b>                   |
| <b>I*</b>    | <b>Filling valve</b>                                    |
| <b>J*</b>    | <b>Filling tube</b>                                     |
| <b>K*</b>    | <b>Seal, PTFE</b>                                       |
| <b>L*/**</b> | <b>SafetyPrime™ valve</b>                               |
| M            | Locking nut                                             |
| <b>N*/**</b> | <b>Backflow tube</b>                                    |
| <b>O*</b>    | <b>Discharge valve (integrated into discharge tube)</b> |
| <b>P*</b>    | <b>Discharge tube</b>                                   |
| <b>Q*</b>    | <b>Screw cap</b>                                        |
| <b>R*</b>    | <b>Seal, PTFE</b>                                       |
| <b>S*</b>    | <b>O-Ring for sterile dispensing (FKM)</b>              |
| T            | Mounting tool                                           |

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## Preparing instrument for use

### Attention!

The adapters supplied with the instrument are made of polypropylene (PP), and can only be applied for media which do not attack PP (► page 36, ETFE-adapters).

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### Warning!


Wear protective clothing, eye protection and gloves. Follow all Safety Instructions and observe the Operating Exclusions and Limitations (► pages 3,4,7,8).

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### Note:

Before using the instrument for the first time, ensure it is rinsed carefully or discard first few samples dispensed.

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## Priming the instrument

### Attention!

To avoid splashes when priming of instruments with backflow valve place toggle in backflow position. When the instrument has no backflow hold discharge tube on inner wall of a suitable recipient to prevent splashes. For small bottles use a bottle stand to prevent tipping over.

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
### Note:


Pull up gently (approx. 30 mm), push down quickly!

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



#### **Warning!**

Always wear protective gloves when touching the instrument and the bottle, especially when using dangerous liquids.  
Carry the mounted instrument only as shown in the figure **2**.

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## Dispensing

### Warning!

Follow all Safety Instructions and observe the Operating Exclusions and Limitations (► pages 4, 8, 9). The discharge tube must always point away from the user or any other person. Never press down the piston when the screw cap is screwed on. Reagent may drip off the discharge tube and screw cap.

### Attention!

After use, always leave the piston in the down position.

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
## Serial dispensing with flexible discharge tube (optional)

### Attention!

Each time you are going to use the tubing, examine it carefully for any damages.  
To dispense aggressive liquids, e.g. concentrated acids and alkalis, you should take safety measures in addition to the normal precautions of wearing protective clothing, eye protection and protective gloves.  
We recommend that you use a protective screen.  
When using the flexible discharge tube, the bottle must be supported using a bottle stand to help prevent the bottle from accidentally tipping over. To help avoid reagent splashing from the tube, always grip the tube firmly by the handle and replace in the holder after use.

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


### **Operating Exclusions:**

The flexible discharge tube is compatible with all liquids suitable for the Dispensette® Organic (except peroxides) and the Dispensette® III (► page 7).

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### **Note:**

The coil of the tubing can be stretched to a maximum length of 800 mm. The entire coil must lie in regular loops and must not be twisted.

When working with the flexible discharge tube, the accuracy and coefficient of variation indicated in the Operating Manual are only obtained for volumes  $\geq 2$  ml and by smooth and steady operation, gently approaching the upper and lower stops. Rinse the tube carefully. Do not dismantle!

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## Connecting drying tube (optional)

### Note:

The screwthreads of the bottle and/or the threaded adapter should be sealed with PTFE tape.

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

## Connecting micro filter connector (optional)

### Note:

The screwthreads of the bottle  
and/or the threaded adapter should  
be sealed with PTFE tape.

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
## Cleaning

The instrument must be cleaned in the following situations to assure correct operation:

- if the piston becomes difficult to move
- after dispensing liquids which form deposits
- before changing the reagent
- prior to long term storage
- should liquid have accumulated in the screw cap
- prior to sterilization
- before replacing the valves

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### Warning!

Follow all Safety Instructions (► page 4-5). The glass cylinder, valves, filling and discharge tubes contain reagent. Point the valves and tube openings away from your body. Wear protective clothing, goggles and appropriate hand protection.

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




### Note:

Never change the pistons of the instruments!

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- 6** Unscrew piston and gently lift it out of the cylinder.  
**A:** for type Digital  
**B:** for type Analog, Fix  
 (➤ page 21)

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## Cleaning

– continued –

- 1** Empty the instrument and unscrew.
- 2** Mount the instrument onto another bottle filled with an appropriate cleaning solution.
- 3** Flush the instrument in dispensing and backflow position.
- 4** Empty the bottle, clean it and fill it with distilled water. Repeat step 3 and 4.
- 5** Empty the instrument.

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**Note:**



Type Digital:  
Assemble and disassemble only at maximum volume.

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- 7 Carefully remove deposits at the edge of the glass cylinder with a screw driver. Clean piston and cylinder with a bottle-brush. Flush instrument with distilled water.

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**Note:**




For type Digital make sure to assemble the dispensing unit in such a way that the catch is inserted below the cylinder!

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


## Autoclaving

This instrument is autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) for 20 minutes in accordance with widely used guidelines. It is the user's responsibility to ensure effective autoclaving.

For frequent autoclaving it is recommended to use filling valves with olive-shaped nozzle (PEEK) to fix the filling tube (► page 38).

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**R** Dichtring, PTFE  
seal, PTFE

**S** O-Ring  
O-ring

## Preparation prior to the first autoclaving

Remove the built-in PTFE-ring (R) and replace with the supplied elastic O-ring (S) made of FKM.

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
### Note:

FKM and PEEK have limited chemical resistance!

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## Autoclaving

Follow all cleaning instructions  
(► page 20).

- 1 + 2** **Loosen** all push-in and screw connections.
- 3** **Tighten** the filling valve.
- 4 – 7** Digital model only: Check that piston mount is firmly seated.
- 8** **Place** all parts (instrument with piston in down position) on a towel and autoclave. Avoid contact with metallic surfaces.

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### Note:

Allow to cool down in the closed autoclave to avoid deformation through too rapid cooling.  
After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them.  
Do not reassemble the instrument until it has cooled down to room temperature.  
(Cooling time approx. 2 hours)

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







**Note:**

Mounting tool can be autoclaved at 121 °C.

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## Checking the volume

The detailed testing instruction can be ordered from BrandTech

- Set to the nominal volume.
- Dispense distilled water.
- Weigh the dispensed quantity on an analytical balance.
- Calculate the dispensed volume taking the temperature into account.
- Perform at least 10 dispensing and weighing operations.
- Calculate the accuracy (A%) and the coefficient of variation (CV%) by means of the formulas used in statistical quality control. The process is described, e.g., in DIN EN ISO 8655-6. Observe the operating manual of the balance manufacturer and the corresponding standards.

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### Calculations (for nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$x_i$  = results of weighings  
 $n$  = number of weighings

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#### Note:

For evaluation and documentation  
 to GLP and ISO we recommend  
 the calibrating software EASYCAL™.

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## Calibration

### Note:

It might be appropriate to calibrate the instrument after a longer period of usage or for specific applications.

1. Calibrate at nominal volume (► page 26).
2. Calculate mean volume (result of weighing) (► page 27).
3. Adjust the instrument (to the calculated mean volume).

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


### Note:

After the adjustment, **no** further calibration check is necessary.

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




### Note:

After the adjustment, it is **necessary** to make another calibration check.

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## Troubleshooting

| Problem                                    | Possible cause                                                                                                                                                                                                                                                                                                                                                                                     | Corrective action                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Piston is difficult to move</b>         | <ul style="list-style-type: none"> <li>- Formation of crystals, dirty</li> </ul>                                                                                                                                                                                                                                                                                                                   | <ul style="list-style-type: none"> <li>- Stop dispensing immediately. Loosen piston by circular motion, but do not disassemble. Follow all cleaning instructions (► page 19).</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Filling not possible</b>                | <ul style="list-style-type: none"> <li>- Volume adjusted to minimum setting</li> <li>- Sticking filling valve</li> </ul>                                                                                                                                                                                                                                                                           | <ul style="list-style-type: none"> <li>- Set to required volume (► page 14).</li> <li>- Clean filling valve (► page 22). If necessary replace the valve with sealing washer.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Instrument does not dispense liquid</b> | <ul style="list-style-type: none"> <li>- Sticking discharge valve</li> </ul>                                                                                                                                                                                                                                                                                                                       | <ul style="list-style-type: none"> <li>- Clean discharge valve. If valve ball is stuck use yellow plastic pipette tip (100 µl) to loosen it. If necessary replace discharge valve.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Air bubbles in the instrument</b>       | <ul style="list-style-type: none"> <li>- Reagent with high vapor pressure has been drawn in too quickly</li> <li>- Seal not inserted, knurled locking nuts not firmly connected</li> <li>- Air not expelled from the instrument</li> <li>- Filling tube not firmly connected or damaged</li> <li>- Valves dirty, not firmly installed or damaged</li> <li>- Backflow tube not connected</li> </ul> | <ul style="list-style-type: none"> <li>- Slowly drag in reagent.</li> <li>- Double check that seal (K) is in place, and that the knurled locking nuts at the discharge tube and the SafetyPrime™ valve are firmly seated and finger tight.</li> <li>- Remove air (► page 13).</li> <li>- Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.</li> <li>- Cleaning procedure (► page 19). Tighten the valves using the mounting tool. If necessary, replace the valves and sealing washers.</li> <li>- Connect backflow tube (► page 10, N).</li> </ul> |
| <b>Dispensed volume is too low</b>         | <ul style="list-style-type: none"> <li>- Discharge tube not firmly connected or damaged</li> <li>- Filling tube not firmly connected or damaged</li> <li>- Valves not firmly connected or damaged</li> </ul>                                                                                                                                                                                       | <ul style="list-style-type: none"> <li>- Push the original manufacturer's discharge tube on firmly. Replace a deformed or damaged discharge tube.</li> <li>- Cleaning procedure (► page 19). Push the filling tube on firmly. If this does not suffice, cut off approx. 1 cm of the tube at the upper end and re-connect it or replace filling tube (► page 37).</li> <li>- Cleaning procedure (► page 19). Tighten the valves using the mounting tool. If necessary, replace valves and sealing washers.</li> </ul>                                                                                                                     |

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**Model Digital • Easy Calibration is  
manufactured under U.S. Patent  
5,957,330.**

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\* Error limits according to DIN EN ISO 8655-5 (A = Accuracy, CV = Coefficient of variation) relative to the nominal capacity (maximum volume) indicated on the instrument, obtained at equal temperature (20 °C) of instrument, ambience and dist. water, and with smooth jerk-free operation.

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## Items supplied

1 Dispensette® III or  
Dispensette® Organic,  
1 performance certificate,  
1 operating manual,  
1 discharge tube  
(with or without SafetyPrime™),  
1 calibrating/mounting tool,  
1 filling tube,  
and adapters of polypropylene for  
following bottle threads:

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## Return for Repair

### Attention!

Transporting of hazardous materials without a permit is a violation of federal law.

BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

- Therefore, contact BrandTech Scientific, Inc. and obtain return authorization **before** sending your instrument for service.

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## Warranty

As provided by law, any and all warranties are null and void if the product has been misused, modified or repaired by unauthorized personnel, if the defects are caused by negligence (instruction manual, maintenance) or by normal wear and tear. Use only original manufacturer's accessory/spare parts.

For conditions and extent of warranty refer to our warranty card enclosed with the product.

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Subject to technical modification without notice. Errors excepted.

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**Note:**

DIN EN ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments.

We recommend to check the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter.

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Neither the Dispensette® III nor the Dispensette® Organic is designed for use with Hydrofluoric acid (HF). Consults the "Operating Exclusions" and "Operating Limitations" sections of the Operating Manual (see pages 7, 8).

The above guidelines have been carefully reviewed prior to publication. Always follow instructions in the operating manual of the instrument as well as reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g. biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BRAND.

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### **Important Product Information**

Filling tubes, discharge tubes and valves are not interchangeable among the several generations of Dispensette® dispensers.

For assistance in selecting the right replacement parts for your newer and older instruments, consult your operating manual, or contact BrandTech Scientific, Inc.

### **Special Upgrade Offer!**

If you have an older Dispensette® dispenser or a dispenser from a different manufacturer, contact BrandTech Scientific for information on how you can upgrade to the Dispensette® III or Dispensette® Organic at special trade-in prices.

For details:

Call BrandTech toll-free: 888-522-2726

Or email: [products@brandtech.com](mailto:products@brandtech.com)