

Discover  
the softer side  
of light pipetting



**AccuPet Pro**

*20 % Softer than Eppendorf Physiocare*

*30 % Lighter than Rainin Pipet-Lite*

*40 % More affordable than all of them †*



# Ergonomically Designed Advance Micropipette

AccuPet Pro™ micropipettes are the new standard in microliter measurements. This line of ergonomically designed continuously adjustable digital pipettes is built with the philosophy of combining attractive form with practical function at an affordable price. The new features allow the user to pipette more comfortably and efficiently yielding results with better accuracy and precision. AccuPet Pro™ micropipettes come in eight standard sizes allowing measurements from 0.2 - 5000 µl. Each pipette includes a GLP/GMP calibration certificate. The user can calibrate pipettes by the simplified EasyCal technique. This highly durable and reliable instruments come with the best pipette guarantee — full three year warranty.

## AccuPet Pro

### Two Models

AccuPet Pro pipettes come in two modifications: Standard blue and white and Color coded for easy volume selection and identification.

### Ergonomic Handle Design


Designed to minimize the amount of hand stress over time lightweight pipettes feature rubberized soft surfaces that help to make pipetting more comfortable than never before. The improved ergonomic handle features a new comfortable grip and “hanger” that allow user to lightly hold pipette during experiments. This minimizes the amount of hand stress over time. Accuracy and precision are improved while reducing the risk of cumulative trauma disorders.

### DirectVue™ Volume Indicator

Monitor volume settings before, during, and after pipetting with the DirectVue™ Volume Indicator. There is no need to use both hands or set samples down to read the volume settings. Volume easily set by turning the plunger without possibility to snag gloves. The user will have the confidence that volume settings are correct and have not drifted over the course of many measurements.

### Universal Tip Fitting Shafts

Fully autoclavable universal shaft is manufactured from the most durable and chemically resistant plastic PVDF which offer highest resistance to shocks and corrosion. It is designed to provide good and consistent seal with most of the standard tips available on the market today. This enables scientists to continue confidently use tips from their favorite supplier.



AccuPet Pro	AccuPet Pro Color		Range (µl)	Accuracy (%)	Precision (%)	List Price (\$)
AP-2	AC-2	Red	0.2-2.0	12.0% - 5.0%	6.0%-1.2%	\$269.00
AP-10	AC-10	Blue	0.5-10.0	5.0% - 1.0%	2.0%-0.8%	\$269.00
AP-20	AC-20	Yellow	2.0-20.0	3.0% - 1.0%	2.0%-0.4%	\$249.00
AP-50	AC-50	Red	5.0-50.0	2.0% - 0.6%	1.0%-0.5%	\$249.00
AP-100	AC-100	Teal	10.0-100.0	3.0% - 0.8%	1.0%-0.2%	\$249.00
AP-200	AC-200	Purple	20.0-200.0	3.0% - 0.8%	1.0%-0.2%	\$249.00
AP-1000	AC-1000	Navy	100.0-1000.0	3.0% - 0.8%	0.6%-0.2%	\$249.00
AP-5000	AC-5000	Teal	500.0-5000.0	1.0% - 0.6%	0.4%-0.2%	\$269.00

† 20% handle surface coverage by soft rubber grip; 30% lighter by weight; 40% lower average retail price.



**Advanced Pipetting Systems**  
 9212 Mira Este Court, Suite 500  
 San Diego • CA 92126  
 800-755-1779 (p) • 858-536-3342 (f)  
 www.pipetpro.com • info@pipetpro.com

Distributed by: