

To order go online or call

**800-243-3232**

**SAVE 30%  
or more**

## Capp *Harmony*

SEROLOGICAL PIPETTES

[cappharmony.capp.dk](http://cappharmony.capp.dk)

- DESIGNED FOR ACCURATE AND RAPID FLUID DISPENSING, THESE POLYSTYRENE PIPETTES ARE STERILE, PYROGEN FREE AND COTTON PLUGGED FOR A VARIETY OF LIQUID HANDLING APPLICATIONS.
- THE SHARP, BLACK, HIGHLY LEGIBLE PERMANENT GRADUATIONS PREVENT AMBIGUITY AND GUARANTEE PRECISE READING AND ACCURATE DISPENSING TO +1-2%. REVERSE GRADUATIONS SHOW VOLUME REMOVED OR REMAINING ON ALL PIPETTES LARGER THAN 1 ML.
- GENEROUS NEGATIVE GRADUATIONS PROVIDE ADDITIONAL PIPETTING VOLUME.
- EACH PIPETTE COMES WITH EFFECTIVE FILTER BARRIER TO PREVENT OVERFLOW OF LIQUID INTO PIPETTE.
- DESIGNED TO FIT ALL AVAILABLE MOTORIZED PIPETTE CONTROLLERS ON THE MARKET.
- PIPETTES ARE SUPPLIED INDIVIDUALLY WRAPPED IN FIBER-FREE COLOR CODED PAPER AND PLASTIC COMBINATION THAT IS EASY TO OPEN. THE UNIFORM MOUTH PIECE FITS ALL MAJOR PIPETTORS AND IS COLOR CODED FOR EASE OF IDENTIFICATION.

**Get Free CappAid  
with purchase of  
5 cases or more**



CAPP HARMONY PIPETTES ARE MADE FROM TRANSPARENT MEDICAL GRADE POLYSTYRENE STERILIZED BY GAMMA IRRADIATION. EACH LOT IS TESTED FOR ENDOTOXINS (C.5 EU/ML) (ISO TO993-11), CYTOTOXICITY (ISO 10993-5) AND ARE GUARANTEED NON-HAEMOLYTIC (ISO 10993-4). IDEAL FOR USE WITHIN TISSUE CULTURE APPLICATIONS.

## Ordering Info

CAT.No.	DESCRIPTION	LIST CASE	SALE CASE	SALE PACK	COLOR
SP-1-C	SEROLOGICAL PIPETTE, 1 ML, YELLOW, 200 PCS/PACK, 800 PCS/CASE	<del>+186</del>	122	36	Yellow
SP-2-C	SEROLOGICAL PIPETTE, 2ML, GREEN, 150 PCS/PACK, 600 PCS/CASE	<del>+178</del>	112	34	Green
SP-5-C	SEROLOGICAL PIPETTE, 5ML, BLUE, 100 PCS/PACK, 400 PCS/CASE	<del>+146</del>	94	26	Blue
SP-10-C	SEROLOGICAL PIPETTE, 10ML, ORANGE, 100 PCS/PACK, 400 PCS/CASE	<del>+156</del>	112	28	Orange
SP-25-C	SEROLOGICAL PIPETTE, 25ML, RED, 50 PCS/PACK, 200 PCS/CASE	<del>+128</del>	98	30	Red
SP-50-C	SEROLOGICAL PIPETTE, 50ML, PURPLE, 45 PCS/PACK, 190 PCS/CASE	<del>+326</del>	196	52	Purple