



**Liquid  
handling  
just got  
easier.**



**costar**



Micrometer head  
 Inner barrel  
 Fixed scale and index line  
 Movable scale  
 Quick setting knob

Dispensing trigger

Tip eject trigger

Body

Ejector bar

Tip holder

Tips

**Related Costar  
products**

**96-Tip rack**

<b>Catalog Number</b>	<b>Volume Range</b>	<b>Color</b>	<b>Quantity/ Case</b>	<b>Sterile</b>
4863	10-200 $\mu$ l	Natural	10 Rack/Case	No
4864	10-200 $\mu$ l	Natural	10 Rack/Case	Yes
4865	10-200 $\mu$ l	Yellow	10 Rack/Case	No

**Bulk replacement tips**

<b>Catalog Number</b>	<b>Volume Range</b>	<b>Color</b>	<b>Quantity/ Case</b>	<b>Sterile</b>
4860	10-200 $\mu$ l	Natural	1000/10 Bags	No
4862	10-200 $\mu$ l	Natural	1000	No
4866	10-200 $\mu$ l	Yellow	1000	No

1 Alewife Center **1-800-492-1110**  
Cambridge, MA  
02140

Printed in USA

**costar**