

WHEATON[®]

CryoELITE[®] Cryogenic Vials



WHEATON CryoELITE[®] Cryogenic Vials offer an unrivaled cap seal for the shipping and storage of precious biosamples. The patented cap seal exceeds DOT and IATA regulations. Optional 2D Data Matrix Bar Code Inserts provide a unique identifier for traceability. The bottom format allows for an unrestricted view of the 2D bar code, convenient for automated scanning. Made from low binding, cryogenic grade virgin polypropylene, with lot certified RNase / DNase and Endotoxin free certification, providing assurance of product integrity. Colored caps allow for color coding in conjunction with WHEATON color coded KeepIT[®] Boxes.

Product Highlights	
<ul style="list-style-type: none">• Unrivaled cap seal exceeds DOT and IATA regulations ensuring ultimate protection of samples during transportation and demanding freeze-thaw handling	<ul style="list-style-type: none">• Lot certified RNase / DNase and Endotoxin Free providing assurance of product integrity
<ul style="list-style-type: none">• Loctagon[™] Vial Skirt provides stability in freestanding position	<ul style="list-style-type: none">• Colored caps allow for color coding projects along with WHEATON colored freezer and storage boxes

CryoELITE® Cryogenic Vials



Freestanding,
Internal Thread



Freestanding,
External Thread















Round Bottom,
Internal Thread



Round Bottom,
External Thread



CryoELITE® Cryogenic Vials, Freestanding

Cat. No.	Size (mL)	Color	Writing Patch	Sterile	Dimensions Dia x H (mm)	Qty / Case
CryoELITE® Cryogenic Vials, Freestanding, Internal Thread						
W985915	1.2	Natural	Yes	Yes	12 x 40	500
W985902	2	Natural	Yes	No	12 x 50	1000
W985903	2	Natural	No	No	12 x 50	1000
W985922	2	Natural	Yes	Yes	12 x 50	500
W985916	2	White	 Yes	Yes	12 x 50	500
W985917	2	Red	 Yes	Yes	12 x 50	500
W985918	2	Pink	 Yes	Yes	12 x 50	500
W985919	2	Yellow	 Yes	Yes	12 x 50	500
W985920	2	Green	 Yes	Yes	12 x 50	500
W985921	2	Blue	 Yes	Yes	12 x 50	500
W985923	3	Natural	Yes	Yes	12 x 63	500
W985924	4	Natural	Yes	Yes	12 x 77	500
W985925	5	Natural	Yes	Yes	12 x 91	500
CryoELITE® Cryogenic Vials, Freestanding, External Thread						
W985862	1.2	Natural	Yes	Yes	12 x 37	500
W985852	2	Natural	Yes	No	12 x 49	1000
W985853	2	Natural	No	No	12 x 49	1000
W985872	2	Natural	Yes	Yes	12 x 49	500
W985863	2	White	 Yes	Yes	12 x 49	500
W985864	2	Red	 Yes	Yes	12 x 49	500
W985865	2	Pink	 Yes	Yes	12 x 49	500
W985866	2	Yellow	 Yes	Yes	12 x 49	500
W985867	2	Green	 Yes	Yes	12 x 49	500
W985868	2	Blue	 Yes	Yes	12 x 49	500
W985869	3	Natural	Yes	Yes	12 x 63	500
W985870	4	Natural	Yes	Yes	12 x 78	500
W985871	5	Natural	Yes	Yes	12 x 93	500

CryoELITE® Cryogenic Vials, Round Bottom

Cat. No.	Size (mL)	Color	Writing Patch	Sterile	Dimensions Dia x H (mm)	Qty / Case
CryoELITE® Vials, Round Bottom, Internal Thread						
W985910	1.2	Natural	Yes	Yes	12 x 39	500
W985911	2	Natural	Yes	Yes	12 x 49	500
W985900	2	Natural	Yes	No	12 x 49	1000
W985901	2	Natural	No	No	12 x 49	1000
W985912	3	Natural	Yes	Yes	12 x 63	500
W985913	4	Natural	Yes	Yes	12 x 76	500
W985914	5	Natural	Yes	Yes	12 x 90	500
CryoELITE® Vials, Round Bottom, External Thread						
W985860	1.2	Natural	Yes	Yes	12 x 35	500
W985861	2	Natural	Yes	Yes	12 x 49	500
W985850	2	Natural	Yes	No	12 x 49	1000
W985851	2	Natural	No	No	12 x 49	1000

For more information on WHEATON products, please visit www.wheaton.com

WHEATON | 1501 North 10th Street | Millville, New Jersey 08332-2038 | USA
800.225.1437 (U.S. & Canada Only) | 856.825.1100 | 856.825.1368 (F)

WSS_027 10/2012